



Product Specification
Slide mailing container, length: 85 mm, cap
attached

Page 1



Product description

Order number	86.972.150
Product description	Slide mailing container, without ventilation, length: 85 mm, cap attached, 50 piece(s)/bag

Product characteristics

Design	without ventilation
Application	Mailing of samples
Example of use	Secondary packaging according to P 650

Size

Length of product	85 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Slide mailing container, length: 85 mm, cap attached

Page 2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	400
Piece(s) / pallet	14400
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0102 kg
Weight of case	4.87 kg
EAN of inner box	4038917123969
EAN case	4038917405928

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT