



**Product Specification**  
**Slide mailing container, length: 85 mm, cap attached**

**Page 1**



**Product description**

Order number	86.972
Product description	Slide mailing container, with ventilation, length: 85 mm, cap attached, 50 piece(s)/bag

**Product characteristics**

Design	with ventilation
Application	Mailing of samples
Example of use	Secondary packaging according to P 650

**Size**

Length of product	85 mm
-------------------	-------

**Material & colours**

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Slide mailing container, length: 85 mm, cap  
attached**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	400
Piece(s) / pallet	14400
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0104 kg
Weight of case	4.97 kg
EAN of inner box	4038917084741
EAN case	4038917405904

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**