



Product Specification
**Polystyrene container, suitable as outer
packaging for cool transport**

Page 1



Product description

Order number	95.1011
Product description	Polystyrene container, suitable as outer packaging for cool transport

Product characteristics

Design	outer packaging for cool transport
Absorbent liner	no

Size

Width of product	71 mm
Height of product	168 mm
Length of product	140 mm

Material & colours

Product material	styrofoam
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Polystyrene container, suitable as outer
packaging for cool transport**

Page 2

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	36
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	36
Piece(s) / pallet	720
Depth of case	584 mm
Width of case	387 mm
Height of case	355 mm
Case volume	0.0802 cbm
Weight of product	0.0371 kg
Weight of case	2.25 kg
EAN of inner box	4038917084987
EAN case	4038917362399

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP
support and is valid without signature.



SARSTEDT