




Product Specification
Polystyrene container, suitable as outer packaging for cool transport

Page 1



Product description

Order number	95.1011
Product description	Polystyrene container, suitable as outer packaging for cool transport

Product characteristics

Design	outer packaging for cool transport
Absorbent liner	no

Size

Width of product	71 mm
Height of product	168 mm
Length of product	140 mm

Material & colours

Product material	styrofoam
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Polystyrene container, suitable as outer
packaging for cool transport**

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 36
Piece(s) / pallet 720
Depth of case 584 mm
Width of case 387 mm
Height of case 360 mm
Case volume 0.0814 cbm
Weight of product 0.0371 kg
Weight of case 2.25 kg
EAN of inner box 4038917084987
EAN case 4038917362399

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT