



Product description

Order number	95.1124
Product description	Cold transport container, french

Product characteristics

Label/ Print	with print
Colour of print/label	red
Absorbent liner	no

Size

Width of product	71 mm
Height of product	170 mm
Length of product	140 mm

Material & colours

Product material	Acrylnitrile butadiene styrene (ABS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	yes

Packaging

Minimum order qty.	200
Piece(s) / inner box	1
Piece(s) / case	6
Piece(s) / pallet	648
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0,0134 cbm
Weight of product	0,3448 kg
Weight of case	2,3 kg
EAN of inner box	4038917085410
EAN case	4038917004091

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

