



### **Product description**

Order number	95.1625
Product description	Cold transport container, flemish

### **Product characteristics**

Label/ Print	with print
Colour of print/label	red
Absorbent liner	no

### **Size**

Width of product	71 mm
Height of product	170 mm
Length of product	140 mm

### **Material & colours**

Product material	Acrylnitrile butadiene styrene (ABS)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	200
Piece(s) / inner box	1
Piece(s) / case	6
Piece(s) / pallet	648
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0,0134 cbm
Weight of product	0,3448 kg
Weight of case	2,3 kg
EAN of inner box	4038917087261
EAN case	4038917005722

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

