



Product description

| | |
|---------------------|-----------------------------------|
| Order number | 95.1625 |
| Product description | Cold transport container, flemish |

Product characteristics

| | |
|-----------------------|------------|
| Label/ Print | with print |
| Colour of print/label | red |
| Absorbent liner | no |

Size

| | |
|-------------------|--------|
| Width of product | 71 mm |
| Height of product | 170 mm |
| Length of product | 140 mm |

Material & colours

| | |
|-------------------|--------------------------------------|
| Product material | Acrylnitrile butadiene styrene (ABS) |
| Colour of product | white |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Minimum order qty. | 200 |
| Piece(s) / inner box | 1 |
| Piece(s) / case | 6 |
| Piece(s) / pallet | 648 |
| Depth of case | 308 mm |
| Width of case | 223 mm |
| Height of case | 195 mm |
| Case volume | 0,0134 cbm |
| Weight of product | 0,3448 kg |
| Weight of case | 2,3 kg |
| EAN of inner box | 4038917087261 |
| EAN case | 4038917005722 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

