



Product description

Order number	95.900
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, flat design, (LxWxH): 198 x 107 x 38 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, Locking mechanism, 50 piece(s)/case
UN identification	yes

Product characteristics

Label/ Print	with print
Colour of print/label	black
Design	flat design
Example of use	external packaging according to P 650/ ADR
Suitable for	diagnostic specimens

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Post transport packaging, 198 x 107 x 38 mm, for
diagnostic specimens**

Page 2

Size

Width of product	107 mm
Height of product	38 mm
Length of product	198 mm

Material & colours

Product material	cardboard
Colour of product	brown

Purity & certification

Satisfies the requirement	P 650, ADR, IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	3000
Depth of case	584 mm
Width of case	366 mm
Height of case	115 mm
Case volume	0.0246 cbm
Weight of product	0.0601 kg
Weight of case	3.44 kg
EAN of inner box	4038917084550
EAN case	4038917064576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP
support and is valid without signature.



SARSTEDT