



Product description

Order number	95.901
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, (LxWxH): 107 x 198 x 50 mm
UN identification	yes

Product characteristics

Label/ Print	with print
Colour of print/label	black
Example of use	external packaging according to P 650/ ADR

Size

Width of product	198 mm
Height of product	50 mm
Length of product	107 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material cardboard
Colour of product brown

Purity & certification

Satisfies the requirement P 650, IATA, ADR
Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 50
Piece(s) / inner box 50
Piece(s) / case 50
Piece(s) / pallet 2100
Depth of case 656 mm
Width of case 414 mm
Height of case 115 mm
Case volume 0,0312 cbm
Weight of product 0,0701 kg
Weight of case 4,04 kg
EAN of inner box 4038917084567
EAN case 4038917064583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

