



Product Specification
**Post transport packaging, 107 x 198 x 50 mm, for
diagnostic specimens**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 95.901 |
| Product description | Post transport packaging, application e.g.: external packaging according to P 650/ ADR, high design, (LxWxH): 107 x 198 x 50 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, Locking mechanism, 50 piece(s)/case |
| UN identification | yes |

Product characteristics

| | |
|-----------------------|--|
| Label/ Print | with print |
| Colour of print/label | black |
| Design | high design |
| Example of use | external packaging according to P 650/ ADR |
| Suitable for | diagnostic specimens |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Width of product | 198 mm |
| Height of product | 50 mm |
| Length of product | 107 mm |

Material & colours

| | |
|-------------------|-----------|
| Product material | cardboard |
| Colour of product | brown |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR, IATA |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 50 |
| Piece(s) / pallet | 2100 |
| Depth of case | 656 mm |
| Width of case | 414 mm |
| Height of case | 115 mm |
| Case volume | 0.0312 cbm |
| Weight of product | 0.0701 kg |
| Weight of case | 4.04 kg |
| EAN of inner box | 4038917084567 |
| EAN case | 4038917064583 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

