



Product description

Order number	95.902
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, flat design, (LxWxH): 220 x 170 x 40 mm, for diagnostic specimens, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, enclosed foam foil pouch, 50 piece(s)/case
UN identification	yes

Product characteristics

Label/ Print	with print
Colour of print/label	black
Design	flat design
Example of use	external packaging according to P 650/ ADR
Suitable for	diagnostic specimens

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Width of product	170 mm
Height of product	40 mm
Length of product	220 mm

Material & colours

Product material	cardboard
Colour of product	brown

Purity & certification

Satisfies the requirement	P 650, IATA, ADR
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	50
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	2400
Depth of case	488 mm
Width of case	401 mm
Height of case	145 mm
Case volume	0,0284 cbm
Weight of product	0,0679 kg
Weight of case	3,86 kg
EAN of inner box	4038917084574
EAN case	4038917069854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

