



Product description

Order number	95.903
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, (LxWxH): 192 x 146 x 75 mm, for cold transport container
UN identification	yes

Product characteristics

Label/ Print	with print
Colour of print/label	black
Example of use	external packaging according to P 650/ ADR
Suitable for	cold transport container

Size

Width of product	146 mm
Height of product	75 mm
Length of product	192 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	cardboard
Colour of product	white

Purity & certification

Satisfies the requirement	P 650, IATA, ADR
Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	yes

Packaging

Minimum order qty.	50
Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	1650
Depth of case	548 mm
Width of case	513 mm
Height of case	115 mm
Case volume	0,0323 cbm
Weight of product	0,0782 kg
Weight of case	4,54 kg
EAN of inner box	4038917084581
EAN case	4038917079594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

