



Product description

Order number	95.997
Product description	Cold transport container for blood gas capillaries, transparent, length: 50 mm

Product characteristics

Suitable for	blood gas capillaries
Absorbent liner	no

Size

Width of product	50 mm
Height of product	140 mm
Length of product	50 mm
Usable length	118 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	transparent
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	66
Piece(s) / pallet	2508
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0,0295 cbm
Weight of product	0,0446 kg
Weight of case	4,59 kg
EAN of inner box	4038917084918
EAN case	4038917362382

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

