



**Product Specification**  
**Mailing containers, length: 104.4 mm, Ø opening:**  
**98 mm, cap assembled**

**Page 1**



**Product description**

Order number	78.9925.417
Product description	Mailing containers, (LxW): 104.4 x 104.4 mm, Ø opening: 98 mm, cap assembled, cap: red

**Product characteristics**

Application	Transportation of diagnostic specimens
Absorbent liner	no

**Size**

Diameter	104.4 mm
Diameter of opening	98 mm
Width of product	104.4 mm
Length of product	104.4 mm
Length including cap	207 mm
Usable length	181 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Mailing containers, length: 104.4 mm, Ø opening:**  
**98 mm, cap assembled**

**Page 2**

**Material & colours**

Product material                      cardboard  
Colour of cap                              red

**Purity & certification**

Latex-free                                 yes  
Pyrogen-free                              no  
Safety data sheet                         no

**Packaging**

Minimum order qty.                      50  
Piece(s) / inner box                      50  
Piece (s) / outer case                    50  
Piece(s) / pallet                         600  
Depth of case                              528 mm  
Width of case                               465 mm  
Height of case                              526 mm  
Case volume                                0.1283 cbm  
Weight of case                              4.8 kg  
EAN of inner box                         4038917155519  
EAN case                                      4038917039420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-02-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**