



**Product Specification**  
**Mailing containers, length: 104.4 mm, Ø opening:**  
**98 mm, cap assembled**

**Page 1**



**Product description**

Order number	78.9925.417
Product description	Mailing containers, (LxW): 104.4 x 104.4 mm, Ø opening: 98 mm, cap assembled, cap: red

**Product characteristics**

Application	Transportation of diagnostic specimens
Absorbent liner	no

**Size**

Diameter	104.4 mm
Diameter of opening	98 mm
Width of product	104.4 mm
Length of product	104.4 mm
Length including cap	207 mm
Usable length	181 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Mailing containers, length: 104.4 mm, Ø opening:  
98 mm, cap assembled**

**Page 2**

**Material & colours**

Product material	cardboard
Colour of cap	red

**Purity & certification**

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

**Packaging**

Minimum order qty.	50
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	600
Depth of case	528 mm
Width of case	465 mm
Height of case	526 mm
Case volume	0.1283 cbm
Weight of case	4.8 kg
EAN of inner box	4038917155519
EAN case	4038917039420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**