



Product Specification
Mailing containers, length: 104.4 mm, Ø opening:
98 mm, cap assembled

Page 1



Product description

Order number 78.9925.417
Product description Mailing containers, (LxW): 104.4 x 104.4 mm, Ø opening: 98 mm, cap assembled, cap: red

Product characteristics

Application Transportation of diagnostic specimens
Absorbent liner no

Size

Diameter 104,4 mm
Diameter of opening 98 mm
Width of product 104,4 mm
Length of product 104,4 mm
Length including cap 207 mm
Usable length 181 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Mailing containers, length: 104.4 mm, Ø opening:
98 mm, cap assembled**

Page 2

Material & colours

Product material cardboard
Colour of cap red

Purity & certification

Latex-free yes
Pyrogen-free no
Safety data sheet no

Packaging

Minimum order qty. 50
Piece(s) / inner box 50
Piece(s) / case 50
Piece(s) / pallet 600
Depth of case 528 mm
Width of case 465 mm
Height of case 526 mm
Case volume 0,1283 cbm
Weight of case 4,8 kg
EAN of inner box 4038917155519
EAN case 4038917039420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT