



Product description

Order number	86.974
Product description	Mailing envelope, material: PE, transparent, (LxW): 440 x 320 mm, 300 piece(s)/case

Product characteristics

Absorbent liner	no
-----------------	----

Size

Width of product	320 mm
Length of product	440 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	300
Type of smallest subpackaging	case
Piece(s) / inner box	300
Piece (s) / outer case	300
Piece(s) / pallet	7200
Depth of case	295 mm
Width of case	320 mm
Height of case	395 mm
Case volume	0.0373 cbm
Weight of product	0.045 kg
Weight of case	14.2 kg
EAN of inner box	4038917084758
EAN case	4038917003537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

