



Product description

Order number	95.1410
Product description	Mailing envelope, material: PE, blue, (LxW): 310 x 180 mm, 50 piece(s)/case

Product characteristics

Absorbent liner	no
-----------------	----

Size

Width of product	180 mm
Length of product	310 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Satisfies the requirement	P 650
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	1000
Piece(s) / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0112 kg
Weight of case	11.5 kg
EAN of inner box	4038917086530
EAN case	4038917005081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

