



Product description

Order number 95.1710
Product description Foam film bag, (LxW): 185 x 230 mm, material: PE

Size

Width of product 230 mm
Length of product 185 mm

Material & colours

Product material Polyethylene (PE)

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 200
Piece(s) / inner box 200
Piece(s) / case 200
Piece(s) / pallet 6400
Depth of case 458 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Foam film bag, (LxW): 185 x 230 mm, PE

Page 2

Width of case	408 mm
Height of case	205 mm
Case volume	0.0383 cbm
Weight of product	0.0018 kg
Weight of case	0.91 kg
EAN of inner box	4038917087490
EAN case	4038917079648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT