



### **Product description**

Order number	95.1713
Product description	Cross-bottom bag, (LxW): 270 x 270 mm, material: PE, 25 piece(s)/inner box

### **Size**

Width of product	270 mm
Length of product	270 mm

### **Material & colours**

Product material	Polyethylene (PE)
------------------	-------------------

### **Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## **Packaging**

Minimum order qty.	25
Type of inner packaging	inner box
Piece(s) / inner box	25
Piece(s) / case	25
Piece(s) / pallet	1800
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0,0126 cbm
Weight of product	0,0582 kg
Weight of case	1,7 kg
EAN of inner box	4038917087513
EAN case	4038917065917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

