



Product description

Order number	95.1713
Product description	Cross-bottom bag, (LxW): 270 x 270 mm, material: PE, for transport case B 17, 25 piece(s)/case

Product characteristics

Suitable for	transport case B 17
--------------	---------------------

Size

Width of product	270 mm
Length of product	270 mm

Material & colours

Product material	Polyethylene (PE)
------------------	-------------------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	25
Type of smallest subpackaging	case
Piece(s) / inner box	25
Piece (s) / outer case	25
Piece(s) / pallet	1800
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0582 kg
Weight of case	1.7 kg
EAN of inner box	4038917087513
EAN case	4038917065917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

