



Product description

Order number 95.1717
Product description Transport case B 17, with lateral rim for block rack,
(LxWxH): 395 x 395 x 435 mm

Product characteristics

Design with lateral rim for block rack

Size

Width of product 395 mm
Height of product 435 mm
Length of product 395 mm

Material & colours

Product material High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	18
Depth of case	403 mm
Width of case	403 mm
Height of case	460 mm
Case volume	0,0747 cbm
Weight of product	3,561 kg
Weight of case	4,5 kg
EAN of inner box	4038917087544
EAN case	4038917065610

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

