



### **Product description**

Order number	95.1720.011
Product description	Transport case I 50, internal case for block rack, (LxWxH): 270 x 170 x 215 mm

### **Product characteristics**

Design	internal case for block rack
Shape of inner tube	square

### **Size**

Width of product	170 mm
Height of product	215 mm
Length of product	270 mm

### **Material & colours**

Product material	Acrylnitrile butadiene styrene (ABS)
------------------	--------------------------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Purity & certification**

Satisfies the requirement	TRBA 100
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	108
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	1.072 kg
Weight of case	1.34 kg
EAN of inner box	4038917159951
EAN case	4038917159944

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

