



## Product description

Order number	95.1720
Product description	Transport case I 100, internal case with attachment, (LxWxH): 360 x 310 x 275 mm

## Product characteristics

Design	internal case with attachment
Shape of inner tube	square

## Size

Width of product	310 mm
Height of product	275 mm
Length of product	360 mm

## Material & colours

Product material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Satisfies the requirement	TRBA 100
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	30
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	1.613 kg
Weight of case	2.05 kg
EAN of inner box	4038917087551
EAN case	4038917076210

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

