



Product description

Order number	95.1725
Product description	Battery bag, without batteries, (LxW): 205 x 115 mm, material: PL, for standard cooling batteries, 12 piece(s)/case

Product characteristics

Design	without batteries
Suitable for	standard cooling batteries

Size

Width of product	115 mm
Length of product	205 mm

Material & colours

Product material	polyester (PL)
------------------	----------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	12
Type of smallest subpackaging	case
Piece(s) / inner box	12
Piece (s) / outer case	12
Piece(s) / pallet	1980
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.4419 kg
Weight of case	5.5 kg
EAN of inner box	4038917087568
EAN case	4038917068376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

