



Product description

Order number	65.676
Product description	Screw cap, natural, suitable for mailing container 85 x 30 mm, 50 piece(s)/case

Product characteristics

Application	Mailing of samples
Suitable for	mailing container 85 x 30 mm

Size

Length of product	15 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap, natural, suitable for mailing container
85 x 30 mm

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 500
Type of inner packaging case
Piece(s) / inner box 50
Piece(s) / case 500
Piece(s) / pallet 39000
Depth of case 418 mm
Width of case 278 mm
Height of case 130 mm
Case volume 0.0151 cbm
Weight of product 0.0028 kg
Weight of case 1.73 kg
EAN of inner box 4038917083300
EAN case 4038917002240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT