



Product Specification
Screw cap, natural, suitable for Mailing container
114 x 44 mm

Page 1



Product description

Order number	65.678
Product description	Screw cap, natural, suitable for Mailing container 114 x 44 mm, 50 piece(s)/case

Product characteristics

Application	Mailing of samples
Suitable for	Mailing container 114 x 44 mm

Size

Length of product	15 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap, natural, suitable for Mailing container
114 x 44 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 250
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 16500
Depth of case 438 mm
Width of case 268 mm
Height of case 155 mm
Case volume 0.0182 cbm
Weight of product 0.0052 kg
Weight of case 1.61 kg
EAN of inner box 4038917083324
EAN case 4038917002264

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

