



Product Specification
Screw cap, natural, suitable for Mailing container
126 x 30 mm

Page 1



Product description

Order number	65.679
Product description	Screw cap, natural, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag

Product characteristics

Application	Mailing of samples
Suitable for	Mailing container 126 x 30 mm

Size

Length of product	20 mm
-------------------	-------

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap, natural, suitable for Mailing container
126 x 30 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 250
Type of smallest subpackaging bag
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 41250
Depth of case 338 mm
Width of case 228 mm
Height of case 120 mm
Case volume 0.0092 cbm
Weight of product 0.0035 kg
Weight of case 1.1 kg
EAN of inner box 4038917083331
EAN case 4038917002271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

