



Product Specification
Mailing container, construction: round, length: 85 mm, Ø opening: 30 mm



Product description

Order number	78.890
Product description	Mailing container, construction: round, length: 85 mm, Ø opening: 30 mm, without cap, 50 piece(s)/inner box

Product characteristics

Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	no
Shape of opening	round
Liquid-tight	yes

Size

Diameter of opening	30 mm
Length of product	85 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Mailing container, construction: round, length: 85 mm, Ø opening: 30 mm

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	natural

Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	500
Type of inner packaging	inner box
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0,0657 cbm
Weight of product	0,0082 kg
Weight of case	4,94 kg
EAN of inner box	4038917084482
EAN case	4038917003292

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-04-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT