



## Product Specification

**Mailing container, transparent, construction:  
round, length: 114 mm, Ø opening: 44 mm,  
without cap**

**Page 1**



### Product description

Order number	78.895
Product description	Mailing container, round, transparent, construction: round, length: 114 mm, Ø opening: 44 mm, without cap, cap: natural, 50 piece(s)/bag

### Product characteristics

Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	no
Shape of opening	round
Liquid-tight	yes

### Size

Diameter of opening	44 mm
Length of product	114 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**

**Mailing container, transparent, construction:  
round, length: 114 mm, Ø opening: 44 mm,  
without cap**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	5000
Depth of case	598 mm
Width of case	398 mm
Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.017 kg
Weight of case	5.25 kg
EAN of inner box	4038917084512
EAN case	4038917003322

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

