



Product description

Order number	77.575.001
Product description	Urine container, 3 l, brown, with light protection, graduated, material: PE, with paper label, english, 30 piece(s)/case

Product characteristics

Label/ Print	with paper label
Graduation	yes
Type of cap	screw cap
Light protection	yes
Application	Urine collection
Viewing window	yes

Size

Sample volume	3.000,00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Urine container, 3 l, brown, with light protection, graduated

Page 2

Length including cap 246 mm
Length excluding cap 243 mm

Material & colours

Product material Polyethylene (PE)
Colour of product brown
Cap material Polypropylene (PP)
Colour of cap yellow

Purity & certification

Product category IVD product
CE certified CE
Safety data sheet nein
Batched yes

Packaging

Minimum order qty. 30
Type of inner packaging case
Piece(s) / inner box 30
Piece(s) / case 30
Piece(s) / pallet 270
Depth of case 643 mm
Width of case 403 mm
Height of case 515 mm
Case volume 0,1335 cbm
Weight of product 0,1136 kg
Weight of case 4,21 kg
EAN of inner box 4038917266222
EAN case 4038917266215

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-22

This document was prepared by EDP support and is valid without signature.



SARSTEDT