



**Product Specification**  
**Urine container, 3 l, with integrated transfer unit,  
brown, with light protection, graduated**

**Page 1**



**Product description**

Order number	77.575.401
Product description	Urine container, Urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with plastic label, english, cap: yellow, with integrated transfer unit, 30 piece(s)/case

**Product characteristics**

Label/ Print	with plastic label
Graduation	yes
Closure type	screw cap
Light protection	yes
Design	with integrated transfer unit
Application	Urine collection
Example of use	closed urine transfer into a V-Monovette® Urine
Volume of work	700 ml
Viewing window	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

### **Material & colours**

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	yellow

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0.1335 cbm
Weight of product	0.1433 kg
Weight of case	5.1 kg
EAN of inner box	4038917319669
EAN case	4038917319652

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

