



Product Specification
**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 77.575.900 |
| Product description | Urine container NFT, Urine collection with hygienically closed, needle-free collection function for Urine Monovettes, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with paper label, German/English, cap: white, with integrated needle-free transfer unit, assembled, 30 piece(s)/case |

Product characteristics

| | |
|-----------------------|---|
| Label/ Print | with paper label |
| Colour of print/label | white/black |
| Graduation | yes |
| Cap | membrane cap |
| Closure type | screw cap |
| Light protection | yes |
| Design | with integrated needle-free transfer unit |
| Application | Urine collection with hygienically closed, needle-free collection function for Urine Monovettes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 2

| | |
|----------------|---|
| Example of use | closed urine transfer into a Monovette® Urine |
| Volume of work | 3.000 ml |
| Viewing window | yes |
| Connection | Luer (m) |

Size

| | |
|----------------------|-------------|
| Sample volume | 3.000.00 ml |
| Max. volume | 3 l |
| Diameter of opening | 83,5 mm |
| Width of product | 130 mm |
| Length of product | 125 mm |
| Length including cap | 246 mm |
| Length excluding cap | 243 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polyethylene (PE) |
| Colour of product | brown |
| Closure material | Polypropylene (PP) |
| Colour of cap | white |

Purity & certification

| | |
|-------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|----------------------------------|--------|
| Minimum order qty. | 30 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 30 |
| Piece (s) / outer case | 30 |
| Piece(s) / pallet | 270 |
| Depth of box | 643 mm |
| Width of box | 403 mm |
| Height of box | 515 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container NFT, 3 l, with integrated
needle-free transfer unit, brown, with light
protection, graduated**

Page 3

| | |
|-------------------|---------------|
| Depth of case | 643 mm |
| Width of case | 525 mm |
| Height of case | 403 mm |
| Case volume | 0.136 cbm |
| Weight of product | 0.1433 kg |
| Weight of case | 5.1 kg |
| EAN of inner box | 4038917947718 |
| EAN case | 4038917947725 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP
support and is valid without signature.



SARSTEDT