



Product description

Order number 77.578
Product description UriSet 24, with stabiliser, 3 l, with transport tube, brown, with light protection, graduated, material: PE, 30 piece(s)/case

Product characteristics

Type of preparation with stabiliser
Graduation yes
Type of cap screw cap
Light protection yes
Design with transport tube
Application Urine collection
Viewing window yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	3.000,00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Cap material	Polypropylene (PP)
Colour of cap	green

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	30
Type of inner packaging	case
Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0,1335 cbm
Weight of product	0,232 kg
Weight of case	7,76 kg
EAN of inner box	4038917082631
EAN case	4038917067751

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

