



## Product Specification

**UriSet 24, with stabiliser, 3 l, with transport tube  
30 ml, brown, with light protection, graduated**

**Page 1**



### Product description

Order number	77.578
Product description	UriSet 24, Complete set for optimum 24-hour urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, with stabiliser, brown, with light protection, graduated, material: PE, cap: green, with transport tube 30 ml, 30 piece(s)/case

### Product characteristics

Type of preparation	with stabiliser
Graduation	yes
Closure type	screw cap
Light protection	yes
Design	with transport tube 30 ml
Application	Complete set for optimum 24-hour urine collection
Viewing window	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**UriSet 24, with stabiliser, 3 l, with transport tube  
30 ml, brown, with light protection, graduated**

**Page 2**

### Size

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

### Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	green

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0.1335 cbm
Weight of product	0.232 kg
Weight of case	7.76 kg
EAN of inner box	4038917082631
EAN case	4038917067751

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**