



Product description

Order number	77.594
Product description	Urine container, 2 l, brown, with light protection, graduated, material: PE, with paper label, french

Product characteristics

Label/ Print	with paper label
Colour of print/label	yellow/black
Graduation	yes
Type of cap	screw cap
Light protection	yes
Application	Urine collection
Viewing window	no

Size

Sample volume	2.000,00 ml
Max. volume	2 l
Diameter of opening	70 mm
Width of product	106 mm
Height of product	240 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Urine container, 2 l, brown, with light protection, graduated

Page 2

Length of product	106 mm
Length including cap	238 mm
Length excluding cap	235 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Cap material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	IVD product
CE certified	CE
Safety data sheet	nein
Batched	yes

Packaging

Minimum order qty.	30
Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	360
Depth of case	563 mm
Width of case	343 mm
Height of case	500 mm
Case volume	0,0966 cbm
Weight of product	0,0956 kg
Weight of case	3,5 kg
EAN of inner box	4038917082761
EAN case	4038917001809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT