



Product description

Order number 84.1402
Product description Label, (LxW): 150 x 55 mm, paper, for urine container 2 and 3 l, german, 30 piece(s)/minigrip bag

Product characteristics

Application Urine collection
Suitable for urine container 2 and 3 l

Size

Width of product 55 mm
Length of product 150 mm

Material & colours

Product material paper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 300
Type of inner packaging minigrip bag
Piece(s) / inner box 300
Piece(s) / case 1800
Piece(s) / pallet 583200
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0,0051 cbm
Weight of product 0,0011 kg
Weight of case 2,28 kg
EAN of inner box 4038917086523
EAN case 4038917078825

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

