



Product description

Order number	75.1354.001
Product description	Container with screw cap, max. working volume: 100 ml, transparent, graduated, material: PP, screw cap, cap: green, cap enclosed, 500 piece(s)/bag

Product characteristics

Graduation	yes
Type of cap	screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter of opening	58 mm
Length excluding cap	73 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	green

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	12000
Case volume	0,0778 cbm
Weight of product	0,0115 kg
Weight of case	7,64 kg
EAN of inner box	4038917247801
EAN case	4038917247795

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

