



Product Specification
Multi-purpose container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent

Page 1



Product description

Order number	75.1355
Product description	Multi-purpose container, max. working volume: 125 ml, (ØxH): 66 x 67 mm, highly transparent, graduated, material: PP, 50 piece(s)/bag

Product characteristics

Graduation	yes
Application	collection and storage of urine
Volume of work	125 ml

Size

Sample volume	125 ml
Max. volume	125 ml
Diameter	66 mm
Diameter of opening	66 mm
Height of product	67 mm
Length excluding cap	67 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Multi-purpose container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	highly transparent

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	24000
Case volume	0,061 cbm
Weight of product	0,0057 kg
Weight of case	6,2 kg
EAN of inner box	4038917409568
EAN case	4038917004879

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT