



### **Product description**

Order number	75.1356
Product description	Urine collection cup, max. working volume: 500 ml, (ØxH): 80 x 148 mm, natural, material: PP, 250 piece(s)/case

### **Product characteristics**

Graduation	no
Application	collection and storage of urine
Volume of work	500 ml

### **Size**

Sample volume	500 ml
Max. volume	500 ml
Diameter	80 mm
Diameter of opening	80 mm
Height of product	148 mm
Length excluding cap	148 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Urine collection cup, 500 ml, (ØxH): 80 x 148 mm,**  
**PP, natural**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product natural

**Purity & certification**

Product category IVD product  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 250  
Type of inner packaging case  
Piece(s) / inner box 250  
Piece(s) / case 250  
Piece(s) / pallet 5000  
Depth of case 458 mm  
Width of case 408 mm  
Height of case 350 mm  
Case volume 0,0654 cbm  
Weight of product 0,03 kg  
Weight of case 8,25 kg  
EAN of inner box 4038917086318  
EAN case 4038917072984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**