



Product description

Order number	75.1356
Product description	Urine collection cup, max. working volume: 500 ml, (ØxH): 80 x 148 mm, natural, material: PP, 250 piece(s)/case

Product characteristics

Graduation	no
Application	collection and storage of urine
Volume of work	500 ml

Size

Sample volume	500 ml
Max. volume	500 ml
Diameter	80 mm
Diameter of opening	80 mm
Height of product	148 mm
Length excluding cap	148 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Urine collection cup, 500 ml, (ØxH): 80 x 148 mm,
PP, natural

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product natural

Purity & certification

Product category IVD product
CE certified CE
Batched yes

Packaging

Minimum order qty. 250
Type of inner packaging case
Piece(s) / inner box 250
Piece(s) / case 250
Piece(s) / pallet 5000
Depth of case 458 mm
Width of case 408 mm
Height of case 350 mm
Case volume 0,0654 cbm
Weight of product 0,03 kg
Weight of case 8,25 kg
EAN of inner box 4038917086318
EAN case 4038917072984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT