



### **Product description**

Order number	75.1437.100
Product description	Urine container, max. working volume: 75 ml, (ØxH): 65 x 48 mm, transparent, graduated, material: PP

### **Product characteristics**

Graduation	yes
Volume of work	75 ml

### **Size**

Sample volume	75 ml
Max. volume	75 ml
Diameter	65 mm
Height of product	48 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Urine container, 75 ml, (ØxH): 65 x 48 mm, PP,**  
**transparent**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 1000  
Piece(s) / inner box 25  
Piece(s) / case 1000  
Piece(s) / pallet 28000  
Case volume 0,065 cbm  
Weight of product 0,0092 kg  
Weight of case 9,65 kg  
EAN of inner box 4038917129138  
EAN case 4038917074711

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**