



Product description

| | |
|---------------------|--|
| Order number | 75.560.901 |
| Product description | Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS |

Product characteristics

| | |
|----------------|--------|
| Graduation | no |
| Volume of work | 148 ml |

Size

| | |
|-------------------|--------|
| Sample volume | 148 ml |
| Diameter | 73 mm |
| Height of product | 67 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | white |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order qty. | 2500 |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 2500 |
| Piece(s) / pallet | 90000 |
| Depth of case | 457 mm |
| Width of case | 368 mm |
| Height of case | 356 mm |
| Case volume | 0.0599 cbm |
| Weight of product | 0.0028 kg |
| Weight of case | 8.16 kg |
| EAN of inner box | 4038917318198 |
| EAN case | 4038917916547 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

