



**Product Specification**  
**Multi-purpose container, 250 ml, (LxØ): 100 x 73 mm, PS, white**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 75.560   |
| Product description | Multi-purpose container, max. working volume: 250 ml, (LxØ): 100 x 73 mm, Ø opening: 73 mm, white, material: PS, 50 piece(s)/bag |

**Product characteristics**

|                |                                 |
|----------------|---------------------------------|
| Graduation     | no                              |
| Application    | collection and storage of urine |
| Volume of work | 250 ml                          |

**Size**

|                      |        |
|----------------------|--------|
| Sample volume        | 250 ml |
| Diameter             | 73 mm  |
| Diameter of opening  | 73 mm  |
| Height of product    | 100 mm |
| Length of product    | 100 mm |
| Length excluding cap | 100 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | white            |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 20000         |
| Depth of case                 | 598 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 165 mm        |
| Case volume                   | 0.0393 cbm    |
| Weight of product             | 0.0076 kg     |
| Weight of case                | 4.3 kg        |
| EAN of inner box              | 4038917339902 |
| EAN case                      | 4038917339896 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

