



Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76
mm, PP, transparent

Page 1



Product description

Order number	75.562.005
Product description	Container with screw cap 100 ml, with graduation, PP, transparent, with assembled lid, length 73 mm, opening diameter 62 mm, 5 pieces in the bag

Product characteristics

Graduation	yes
Type of cap	Screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Container with screw cap, 100 ml, (ØxH): 57 x 76
mm, PP, transparent

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	white

Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	200
Type of inner packaging	bag
Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0,0597 cbm
Weight of product	0,0153 kg
Weight of case	3,74 kg
EAN of inner box	4038917794053
EAN case	4038917794060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-10-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT