



## Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown

Page 1



### Product description

Order number

75.562.011

Product description

Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, brown, with light protection, graduated, material: PP, screw cap, cap: brown, cap enclosed, 50 piece(s)/bag

### Product characteristics

Graduation

yes

Closure type

screw cap

Light protection

yes

Application

collection and storage of urine

Volume of work

100 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown**

**Page 2**

**Size**

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	7500
Depth of case	443 mm
Width of case	323 mm
Height of case	330 mm
Case volume	0.0472 cbm
Weight of product	0.0153 kg
Weight of case	4.54 kg
EAN of inner box	4038917113618
EAN case	4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

