



### Product description

Order number	75.562.011
Product description	Container with screw cap, max. working volume: 100 ml, (ØxH): 57 x 76 mm, brown, with light protection, graduated, material: PP, screw cap, cap: brown, cap enclosed, 50 piece(s)/bag

### Product characteristics

Graduation	yes
Type of cap	screw cap
Light protection	yes
Application	collection and storage of urine
Volume of work	100 ml

### Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown**

**Page 2**

Height of product 76 mm  
Length excluding cap 72 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product brown  
Cap material High Density Polyethylene (HD-PE)  
Colour of cap brown

**Purity & certification**

Product category IVD product  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 250  
Type of inner packaging bag  
Piece(s) / inner box 50  
Piece(s) / case 250  
Piece(s) / pallet 7500  
Depth of case 443 mm  
Width of case 323 mm  
Height of case 330 mm  
Case volume 0,0472 cbm  
Weight of product 0,0153 kg  
Weight of case 4,54 kg  
EAN of inner box 4038917113618  
EAN case 4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**