



## Product Specification

Urine container NFT, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated needle-free transfer unit, transparent

Page 1



### Product description

Order number	75.562.900
Product description	Urine container NFT, Urine collection with hygienically closed, needle-free collection function for Urine Monovettes, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, with safety label, with integrated needle-free transfer unit, cap assembled, sterile, 5 piece(s)/bag

### Product characteristics

Type of collection	Aspiration
Label/ Print	with safety label
Colour of print/label	white/black
Graduation	yes
Cap	screw cap with needle-free transfer device
Closure type	screw cap
Light protection	no
Design	with integrated needle-free transfer unit

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine container NFT, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated needle-free transfer unit, transparent**

**Page 2**

Application	Urine collection with hygienically closed, needle-free collection function for Urine Monovettes
Example of use	closed urine transfer into a Monovette® Urine
Suitable for	Urine Monovette®
Volume of work	100 ml
Connection	Luer (m)

### Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	200
Piece(s) / pallet	5600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine container NFT, 100 ml, (ØxH): 57 x 76 mm,  
PP, with safety label, with integrated needle-free  
transfer unit, transparent**

**Page 3**

Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0186 kg
Weight of case	4.58 kg
EAN of inner box	4038917934176
EAN case	4038917934183

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**