



Product description

| | |
|---------------------|--|
| Order number | 75.563 |
| Product description | Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø: 57 mm, Ø opening: 57 mm, transparent, graduated, material: PP, screw cap, without cap, 50 piece(s)/bag |

Product characteristics

| | |
|----------------|---------------------------------|
| Graduation | yes |
| Closure type | screw cap |
| Application | collection and storage of urine |
| Volume of work | 100 ml |

Size

| | |
|----------------------|--------|
| Sample volume | 100 ml |
| Max. volume | 100 ml |
| Diameter | 57 mm |
| Diameter of opening | 57 mm |
| Length of product | 72 mm |
| Length excluding cap | 72 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 15000 |
| Depth of case | 443 mm |
| Width of case | 323 mm |
| Height of case | 330 mm |
| Case volume | 0.0472 cbm |
| Weight of product | 0.0084 kg |
| Weight of case | 5.01 kg |
| EAN of inner box | 4038917082501 |
| EAN case | 4038917001588 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

