



Product description

Order number	75.563
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø: 57 mm, Ø opening: 57 mm, transparent, graduated, material: PP, screw cap, without cap, 50 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	57 mm
Length of product	72 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	15000
Depth of case	443 mm
Width of case	323 mm
Height of case	330 mm
Case volume	0.0472 cbm
Weight of product	0.0084 kg
Weight of case	5.01 kg
EAN of inner box	4038917082501
EAN case	4038917001588

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

